

Tune to your favorite channel...



- Ion Channel Cell Lines & Assays - Purinergic Receptor P2X Ion Channels Panel

We offer a large P2X Ion Channels Panel for your Hit Discovery programs

Target Name	Species	Fluorescence Assay	Luminescence Assay	Manual Patch Clamp	Automatic Patch Clamp
P2RX1	Human	*	✓	✓	*
P2rx1	Mouse	*	✓	*	*
P2RX2/P2RX3	Human	✓	*	*	*
P2RX3	Human	*	✓	*	✓
P2rx3	Rat	*	✓	*	*
P2RX4	Human	✓	*	*	✓
P2rx4	Mouse	✓	*	*	✓
P2RX7	Human	✓	*	*	✓
P2rx7	Mouse	✓	*	*	✓
P2rx7	Rat	✓	*	*	✓

✓ Available

* Please inquire for availability

(release November 2017)

- We can provide you with off-the shelf validated functional cell lines & assays
- The assays are fully validated for high throughput screening (HTS) and compound profiling activities
- We can screen your target of interest against a library provided by you or against our AXX^{DIV3.0} library of 275,000 compounds (or against a subset of the same)

Contact us to receive further information on our products and services

Axxam SpA
Open Zone
Via Meucci 3 – 20091, Bresso, Milan (Italy)

Phone +39 02 2105634
info@axxam.com
www.axxam.com

AXXAM
INNOVATES, ENABLES, PERFORMS

Please be aware that the provision of the above listed assays may require third party licenses. All products are for research use only.